

Thanks for choosing an Activex Weighing Scale

1. Calibrating your scale:

- 1.1 Press the platform centre to turn ON the scale and remove your foot.
- 1.2 "0.0" will be displayed. This position is now calibrated as the 0.0 of the scale.
- 1.3 The scale will switch off and is now ready for use.

** Remember that you'll have to do this every time you lift the scale off the ground and keep it somewhere else. Failing this you will see a 2Kg overweight on the reading because of the scale's own weight.*

2. Using the Scale

- 2.1 Step on and stand very still while the scale computes your weight.
- 2.2 Your weight is displayed and fixed after flashing for three times.
- 2.3 Step off. Your weight will be displayed for a few seconds.
- 2.4 The scale will switch off Automatically.

3. Indicator:

"Err": The accurate working range of the scale is from 10 Kgs – 180 Kgs. If a weight of over 182 kgs is placed on the scale, it would show an Error message **"Err"**. Please get off the scale to avoid causing permanent damage to the load cell.

"Lo": Low Power indicator, If the scale shows **"Lo"** it's time to change the batteries, the scale takes 1 CR2032 battery.

4. Warranty

The scale has a **6 month warranty** out of the box from the date of purchase. We will gladly double this warranty for 2 minutes of your time.

Fill the survey below to extend your warranty to 1 Year.

www.MyActivex.com/survey

For full text of our warranty text
and terms go to
www.MyActiveX.com/terms

Name of Importer: IOSX Trading Pvt. Ltd
Address of Importer: E 1003, Nancy Lake,
Opp. Bharti Vidyapeeth,
Satara Road, Katraj, Pune- 411046
Email ID: support@myactivex.com
Website: MyActiveX.com

ANY QUESTIONS ? Write to Us...
Support@MyActiveX.com





ACTIVE X